

Bacteria Genomic DNA Extraction Kit

Instruction for Use (V1.1)

[REF] UBE-S00350Q

[Specification] 50 preps/kit

[Research Use Only]

Ultrassay Bacteria Genomic DNA Kit uses lysozyme and moderate lysis buffer to lyse cells. Proteinase K is used for protein digestion and RNase A used for RNA digestion. DNA is specifically bound to silica-based column in hypersaline condition, and DNA is eluted by low salt and high pH solution. This kit is suitable for isolating high quality genomic DNA from Gram-positive and Gram-negative bacteria. The isolated DNA is suitable for PCR, restriction enzyme digestion, and Southern blot.

- Fast: the whole process can be completed in 50 minutes
- High yield: DNA yield up to 20 µg

[Kit Contents]

No.	Content	50 preps
1	Resuspension Buffer 11 (RB11)	12ml
2	Lysis Buffer11 (LB11)	6ml
3	Binding Buffer11 (BB11)	10 ml
4	Clean Buffer 11 (CB11)	55 ml
5	Wash Buffer 11 (WB11)	12 ml
6	Elution Buffer (EB)	25 ml
7	RNase A (10 mg/ml)	1 ml
8	Proteinase K (20 mg/ml)	1 ml
9	Genomic Spin Columns with Collection Tubes	50 pcs

[Storage Conditions and Shelf Life]

RNase A at -20°C for 24 months, others at room temperature (15-25°C) for 12 months.

[Acceptable Samples]

Gram-positive or Gram-negative bacteria cells ≤109.

[Test Procedures]

Before use, add appropriate volume of 96-100% ethanol to BB11 (15ml) and WB11 (48ml).

- Lysozyme will be supplied by users. Prepare fresh lysozyme/RB11 mix for each use (4 mg lysozyme/ 200 µl RB11)
- Prepare 70% ethanol for extraction of Gram-positive coccus; prepare glass bead for breaking Actinomyces hyphae clump.

(All centrifugation steps are performed at room temperature)

1. Material treatment

Lysis of Gram-negative Bacteria

- Transfer 1 ml of overnight Gram-negative bacteria to a 1.5 ml tube and centrifuge at 12,000×g for 1 minute. Discard the supernatant.
- Add 100 µl of LB11 and 20 µl of Proteinase K into the tube. Resuspend the bacteria by vortex.
- Incubate at 55°C for 15 minutes. (Solution should be clear after incubation. If not, extend the incubation time to 30 minutes, vortex for every 5 minutes.)










Lysis of Gram-positive Bacteria

- Transfer 1 ml of overnight Gram-positive bacteria to a 1.5 ml tube and centrifuge at 12,000×g for 1 minute. Discard the supernatant.
(Note: when extract Gram-positive coccus, resuspend it with 500 µl of 70% ethanol, incubate on ice for 20 minutes and then centrifuge at 10,000×g for 1 minute, discard the supernatant, then process to step (b). When extract Actinomyces, use glass bead to break the hyphae clump, centrifuge at 10,000×g for 1 minute, discard the supernatant, then process to step (b).)
 - Resuspend the bacteria by adding 200 µl of RB11 (containing 4 mg lysozyme) to the tube. Incubate at 37°C for at least 60 minutes (Note: the incubation time can be extended to 3 hours if large amount of bacteria is used), and centrifuge at 10,000×g for 1 minute. Discard the supernatant.
 - Add 100 µl of LB11 and 20 µl of Proteinase K into the tube. Resuspend the bacteria by vortex.
 - Incubate at 55°C for 15 minutes. (Solution should be clear after incubation. If not, extend the incubation time to 30 minutes, vortex for every 5 minutes.)
- Add 20 µl of RNase A to the tube, mix and incubate at room temperature for 2 minutes.
 - Add 400 µl of BB11 (check to make sure 96-100% ethanol has been added) and vortex for 30 seconds. (White flocculent precipitate or transparent gelatinous matter may present in this step, this would not affect DNA extraction)
 - Transfer the entire contents to a spin column, centrifuge at 12,000×g for 30 seconds, discard the flow-through.
 - Add 500 µl of CB11, centrifuge at 12,000×g for 30 seconds, and discard the flow-through.
 - Repeat step 5 once.
 - Add 500 µl of WB11 (check to make sure 96-100% ethanol has been added), centrifuge at 12,000×g for 30 seconds, discard the flow-through.
 - Repeat step 7 once.
 - Centrifuge at 12,000×g for 2 minutes to remove residual WB11.
 - Place the spin column in a sterile 1.5 ml microcentrifuge tube. Add 50-200 µl of Elution Buffer (preheated to 65°C) or sterile, distilled water (pH >7.0) to the center of column. Incubate at room temperature for 2 minutes. Centrifuge at 12,000×g for 1 Lysis of Gram-positive Bacteria minute to elute genomic DNA.
 - Repeat step 10 once. Store the isolated DNA at -20°C.

Notes

- To avoid incomplete lysis, do not use too much starting materials.
- Use sterile tubes and pipette tips to avoid contaminations.

[Index of Symbols]

Symbols	Meanings	Symbols	Meanings
	Use By		Date of Manufacture
	Temperature Limitation		Consult Instructions for Use
	Lot Number		Manufacturer
	Number of Tests		Reference Number
			Any warnings and/or precaution to take



Ultrassay Biotech Co., Ltd.

Add: 2906, Building B, Dongyi Square, No. 169 Funan Road, Hefei, Anhui, China

Tel: +86-551-6288 1663

Website: [Http://ultrassay.com](http://ultrassay.com)

Version: A/0 Issue date: 2021.12.21