



## Pipetting function

Parameter Knob press → Pipetting,  
Long press → Mixing




## Stepper function

Maximum steps = Pipette nominal volume /  
One step dispensing volume



## Other function

Pipetting speed adjustment  
Key tone ON/OFF



## Supports dual charging mode

### USB Charging



### Support automatic Calibration

Connect it to the computer, with the free software  
provided for calibration

## Contact Charging



Complimentary Pipette Holder

# Technical Specification

8-channel Adjustable Volume Electronic Pipettes						
Volume Range	Increment	Test Volume	Systematic Error		Radomn Error	
µL	µL	µL	µL	%	µL	%
0.5~10	0.01	10	±0.2	±2.0	±0.10	±1.0
		5	±0.2	±4.0	±0.10	±2.0
		1	±0.08	±8.0	±0.05	±5.0
10~100	0.1	100	±0.8	±0.8	±0.3	±0.3
		50	±0.5	±1.0	±0.4	±0.8
		10	±0.3	±3.0	±0.2	±2.0
15~300	1	300	±1.8	±0.6	±0.9	±0.3
		150	±1.5	±1.0	±0.75	±0.5
		15	±0.6	±4.0	±0.015	±1.0
30~300	1	300	±1.8	±0.6	±0.9	±0.3
		150	±1.5	±1.0	±0.75	±0.5
		30	±0.9	±3.0	±0.63	±1.0

\*Ultrassay BioTech Co., Ltd. reserves the right of final interpretation of different application.

# dPette+

## Multi-functional 8-channel Electronic Pipette

Multi-functional high performance 8-channel electronic pipette offers productive pipetting with its easy to understand operations. Its uniquely lightweight streamlined design ensures an effortless transfer of multiple samples with an increased throughput and data reproducibility.



### • Easy Operation

Intuitive interface for setting functions and parameters.

### • Ergonomic Design

Low operation forces for complete work, which is exceptionally fatigue-free.

### • Convenient and Versatile

360° pipetting. Double knobs for simple and versatile control. Easy loading tip cones offers smooth and leak free operation.

### • Dual charging modes

Use of USB charger or the charging stand to ensure uninterrupted use.

## Three easy steps to operate dPette+

1

Long press the Parameter Knob for 2 seconds to Start.

2

Rotate the knob back and forth quickly to Switch Pipetting, Continuous Dispensing, and other Function Settings

3

Quickly turn the Parameter Knob to unlock it, turn it to adjust the parameter such as volume, press it to pipette, long press it to mix

