



## Mini Vortex Mixer

### Introduction

Ultrassay mini vortex mixer is an ideal vortex mixer with its tiny size and easy-to-use, and it is designed for oscillating and blending of various reagents, solutions, chemicals and other samples. Simply push the tube into the soft attachment to start shaking. It is used for mixing samples in small containers (0.2 ~ 50mL and diameter less than 108mm tubes) and has the function of inching and continuous.



### Feature

- Two operating modes: Inching & Continuous;
- Exquisite appearance: Tiny size, advanced design, freely use;
- Corrosion-resistant material: Adopt to special rubber, strong corrosion-resistance;
- Supporting outdoor: Mobile power pack is optional.

Product Parameters		
Model	UVM-Mini	UVM-Mini+
Operation & Display	/	/
Mixing Motion	Orbital	Orbital
Orbital Diameter	4.5mm	4.5mm
Operation Mode	Inching	Inching & Rotary button
Rotational Speed	0 ~ 4200rpm	0 ~ 4200rpm
Timing Range	/	/
Power Supply	220V/50Hz	220V/50Hz
Power Consumption	7W	7W
External Dimension	W102 x D102 x H88mm	W102 x D102 x H88mm
Net Weight	0.8kg	0.8kg