

# DISPOSABLE CERVICAL SAMPLER



## Description

Disposable cervical sampler (Type A/B) is used for clinical collection of cervical cells for TCT liquid based cytology tests etc., It is designed according to the anatomical structure of women's cervix, the shape fits the cervical anatomy. The product is composed of PE brush head, PP rod and pushing tube, which is used for clinical sampling of gynecological cervical examination, With soft texture, no damage to the cervical mucosa. Using the push tube design, the brush head is pushed into the preservation medium bottle after sampling to avoid secondary contamination.

## Order Information

| Item No.  | Product Name                | Specification | Packing Size            |
|-----------|-----------------------------|---------------|-------------------------|
| UB-SC0100 | Disposable Cervical Sampler | Type A        | 50pcs/box, 1000pcs/case |
| UB-SC0101 | Disposable Cervical Sampler | Type B        | 50pcs/box, 1000pcs/case |

**Ultrassay BioTech Co., Ltd.**

Add: RM3018 Poly D9, No. 999 Luzhou Avenue, Baohe District, Hefei, Anhui, China

Tel: +86-551-62881663

E-mail: [info@ultrassay.com](mailto:info@ultrassay.com)

Website: <https://ultrassay.com>



# DISPOSABLE CERVICAL SAMPLER



## Description

Disposable cervical sampler (Type C) is used for HPV examination, With soft texture, no damage to the cervical mucosa. It can collect more tissue samples, and simple operation and analysis.

The breakpoint is designed at the rod of the type C sampler. And after sampling, the rod is broken to facilitate the brush head to be stored in the preservation medium for further process. Blister packing, EO sterile.

## Order Information

| Item No.  | Product Name                | Specification             | Packing Size            |
|-----------|-----------------------------|---------------------------|-------------------------|
| UB-SC0102 | Disposable Cervical Sampler | Type C (5.5cm breakpoint) | 50pcs/box, 1000pcs/case |
| UB-SC0103 | Disposable Cervical Sampler | Type C (7.0cm breakpoint) | 50pcs/box, 1000pcs/case |

**Ultrassay BioTech Co., Ltd.**

Add: RM3018 Poly D9, No. 999 Luzhou Avenue, Baohe District, Hefei, Anhui, China

Tel: +86-551-62881663

E-mail: [info@ultrassay.com](mailto:info@ultrassay.com)

Website: <https://ultrassay.com>



# DISPOSABLE CERVICAL SAMPLER



UB-SC0118



UB-SC0119



UB-SC0120

## Description

The disposable cervical cell sampler (self-sampling kit) is used for women to collect cervical cell samples at home.

The disposable cervical cell sampler (self-sampling kit) is composed of protected shell, rotatory rod, brush head and cap head.

The samplers are divided into A, B and C types according to the style of brush heads. Type A is made of medical grade nylon flocked, Type B is made of medical grade PE, Type C is made of medical grade sponge. EO sterile.

## Order Information

| Item No.  | Product Name                | Specification | Packing Size |
|-----------|-----------------------------|---------------|--------------|
| UB-SC0118 | Disposable Cervical Sampler | Type A        | 800pcs/case  |
| UB-SC0119 | Disposable Cervical Sampler | Type B        | 800pcs/case  |
| UB-SC0120 | Disposable Cervical Sampler | Type C        | 800pcs/case  |

**Ultrassay BioTech Co., Ltd.**

Add: RM3018 Poly D9, No. 999 Luzhou Avenue, Baohe District, Hefei, Anhui, China

Tel: +86-551-62881663

E-mail: [info@ultrassay.com](mailto:info@ultrassay.com)

Website: <https://ultrassay.com>





# HPV SAMPLING KIT



## Description

The HPV Sampling Kit consists of a disposable cervical sampler Type C and a preservation tube (containing the preservation medium).

The HPV-2ML contains a 5ml preservation tube with 2ml preservation medium;

The HPV-3ML contains a 7ml preservation tube with 3ml preservation medium.

## Order Information

| Item No.  | Product Name     | Specification | Packing Size           |
|-----------|------------------|---------------|------------------------|
| UB-SC0115 | HPV Sampling Kit | HPV-2ml       | 50pcs/box, 500pcs/case |
| UB-SC0116 | HPV Sampling Kit | HPV-3ml       | 50pcs/box, 500pcs/case |

**Ultrassay BioTech Co., Ltd.**

Add: RM3018 Poly D9, No. 999 Luzhou Avenue, Baohe District, Hefei, Anhui, China

Tel: +86-551-62881663

E-mail: [info@ultrassay.com](mailto:info@ultrassay.com)

Website: <https://ultrassay.com>

