

dPette+

Multi functional Electronic Pipette

dPette

Simple Electronic Pipette

- Motor driven digitally control pipette with multifunctions
- Easy Operations and Handling
- Intuitive menu interface settings: functions and operations
- Minimal force for pipetting operation
- High performance ensuring accuracy & repeatability
- 2 buttons for all operational settings
- Adjustable speed for aspiration and dispensing
- Li-ion battery enable longer operation time
- Convenient and Versatile
- Self-calibration applicable



Technical Specification

Single-channel Adjustable Volume Electronic Pipettes						
Volume Range	Increment	Test Volume	Systematic Error		Radomn Error	
μL	μL	μL	μL	%	μL	%
0.5~10	0.01	10	± 0.1	± 1.0	± 0.05	± 0.5
		5	± 0.1	± 2.0	± 0.10	± 2.0
		1	± 0.05	± 5.0	± 0.03	± 3.0
5~50	0.1	50	± 0.4	± 0.8	± 0.15	± 0.3
		25	± 0.25	± 1.0	± 0.25	± 1.0
		5	± 0.2	± 4.0	± 0.125	± 2.5
30~300	1	300	± 1.8	± 0.6	± 0.6	± 0.2
		100	± 1.0	± 1.0	± 0.4	± 0.4
		30	± 1.2	± 4.0	± 0.21	± 0.7
100~1000	5	1000	± 6.0	± 0.6	± 2.0	± 0.2
		500	± 5.0	± 1.0	± 1.0	± 0.2
		100	± 3.0	± 3.0	± 0.6	± 0.6

*Ultrassay BioTech Co., Ltd. reserves the right of final interpretation of different application.

dPette/dPette+

Electronic Pipette

dPette is an innovative Electronic Pipette developed. It combines manual pipette features, as ergonomics and lightweight, with electronic pipette features, such as labor-saving and high accuracy, offering new pipetting experience to users.

